



DILUWEL®

Increase your productivity  
quality  
traceability

### PERFORMANCES



Weighing :  
1 g to 4000 g



225ml distribution :  
- 12 seconds  
- 9 seconds (turbo)



> 98 %

# DILUWEL®



**Finger-stop nozzle**  
Avoids cross contamination  
during installation



**Fast & accurate!**  
225ml dispensed  
in 12 seconds (9 in turbo)  
with >98-99% accuracy !



**Safety cabinet use**



**PATENTED** vertical weighing tray  
(3 points for higher stability and  
higher weighing accuracy)

Removable 1 to 4 pumps  
Plug & play - Flexible  
Free positioning = save space



**Compact and light**  
8 kg  
50% more compact  
than competition models!



0°  
Built-in bubble level  
and adjustable legs



**eTRACE®**  
OPEN GATE TO ALL CONNECTIONS



«LABPAGE»  
Full LIVE data tracing



**Watertight to spills**  
Perfect cleanability



CONNECTED  
MICROBIOLOGY





## SPECIFICATIONS

# DILUWEL®



### PERFORMANCES

Fast & accurate dispensing	12 sec. for 225ml / (9 sec in turbo mode)
Dilution accuracy	> 98%
Turbo mode	✓
Bubble level integrated	✓
Motorized dispensing nozzle	✓
Editable programs	19
Weighing range:	1 g to 4000 g
Dilution factor	1.00 to 1/100.00
Anticipation mode	✓
Total dosage mode	✓

### BENEFITS

Stainless steel bag holder – extended lifetime	✓
Watertight base – easy to clean	✓
Motorized dispensing nozzle	✓
19 editable programs	✓
Data management (diluent, users)	✓
2 USB port: 1 for upgrades / 1 for communication module	✓
“Fingerstop” large dispensing nozzle	✓
Plug and play removable pumps	✓
Interactive display (color codes)	✓

### DIMENSIONS & ELECTRICAL SUPPLY

Compacity W x D x H	250 x 275 x 400 mm
Weight	8 kg
Electrical supply / Power	230 V – 50Hz / 120W