



**Retire your
test tubes.**

New: Inlabtec Serial Dilution System.


Faster and up to ten times less expensive than serial dilutions using test tubes.

Automated **Serial Diluter** dispenses diluent accurately and reliably to the sample.

High purity disposable **Serial Dilution Bags** eliminate the laborious preparation of test tubes.

Simply by pressing a button, samples in the **Serial Dilution Bags** are rapidly and homogeneously diluted.

At the end of the process - used bags are simply placed in the **Bag Shell** or can be directly disposed of.



The Serial Diluter is an innovative device to enable simple, fast and accurate execution of serial dilutions for bacterial counts in food safety testing.

**Order a test device now:
www.inlabtec.com**

Saves you time by eliminating

- washing, filling, sealing and autoclaving of test tubes
- need for volume control after autoclaving
- need for storage of the prepared dilution tubes
- manual handling of test tubes and caps during pipetting
- manual mixing

Key benefits are

- Great time savings by eliminating the preparation of test tubes
- Improved consistency by automatic mixing and dilution
- Enhanced workflow efficiency from ergonomic product design and methodology
- Cost savings even at low sample numbers
- Short payback time on initial investment

Distributor:

 **INLABTEC**

INLABTEC AG Oberstraße 149 CH-9000 St. Gallen
+41 71 222 48 65 info@inlabtec.com www.inlabtec.com